

SecurePortIV[®]

Catheter Securement Adhesive
By H.B. Fuller Medical Adhesive Technologies

Introduction to SecurePortIV[®] Evidence-Based Strategies for Preventing Vascular-Access Infections and Maximizing Value

Jennifer Perou, RN
Regional Business Manager-CA/HI
562.619.3022
Jennifer.Perou@hbfuller.com



What is SecurePortIV[®]

- Cyanoacrylate Tissue adhesive that **secures the catheter** and seals the insertion site to **reduce bleeding and oozing** and to **protect the insertion site from microbial contamination**.

How do we support it:

- Education and in-service provided live and/or virtually by a team of **expert Vascular Access Specialists and Professionals** with decades of experience.

SecurePortIV[®] Catheter Securement Adhesive:

Class II Device, 510(k) # K170505


- **First and only** FDA-cleared tissue adhesive with IV catheter use indications
- **INS 2021 and 2024 Infusion Therapy Standards of Practice** with multiple practice recommendations for securement, hemostasis and infection prevention.


Indications for Use

- SecurePortIV® to be applied as a film forming securement and sealant at the point of vascular access catheter skin entry.
- protect the catheter skin entry site by creating a sealant that
 - **immobilizes** surface **bacteria**, preventing them from entering the catheter skin entry site
 - provides a **moisture barrier**
 - used with a transparent film dressing
 - for the **securement** of short-term and long-term vascular access catheters

Contraindications

- Do not apply on wounds with evidence of microbial, bacterial, or fungal infections or gangrene
- Do not apply on mucosal surfaces, dense hair, or skin regularly exposed to body fluid
- Do not use on patients with hypersensitivity to cyanoacrylates or formaldehyde.





VASCULAR ACCESS CHALLENGES



**Bleeding and
Oozing**



**Migration
Dislodgement**



**High failure
rate**




Infection




**Dressing failure
Early dressing
disruption**

Value Proposition



2-octyl cyanoacrylate (80%) butyl cyanoacrylate (CA)



Water diffusion from microbe cell into CA causes cell lysis and death.

Low water activity is a limiting factor for microbial survival

Completely kills ≥ 8 log of bacteria after a 3-minute contact time

Broad spectrum:
Gram –, Gram +, Yeast, Fungi; 3.47-log *C. auris*

SECUREMENT Standard 36

Tissue Adhesive


Medical-grade
cyanoacrylate (CA) glue

SecurePortIV®

Catheter Securement Adhesive
Highly Purified Medical Cyanoacrylate

Subcutaneous Anchor

Securement System
(SASS)




Integrated Securement Device (ISD)

A device that combines
a dressing with
securement functions

Adhesive Securement Device (ASD)

An adhesive-backed
device that adheres to
the skin

INS Standards of Practice



- **Securement technology** for all IV catheter types - **Level I**
- Provides **enhanced catheter stabilization** for PIVCs - **Level II**
- **Reduces PIVC complication rates** - **Level II**
- Use as an alternative to sutures to **reduce** complications, including **infection and dislodgement** - **Level I**
- **Reduces failure and improves dwell time** when used with sutures - **Level III**
- Use in addition to primary dressing **reduces PIVC failures** - **Level III**
- **Promotes decreased dressing changes and increased dwell time** for all IV catheter types - **Level III**
- **Improves hemostasis** and **reduces early dressing changes** - **Level II**
- **Consider use in the neonatal population** - **Level V**
- **Reduces microorganism growth** on the catheter tip - **Level IV**
- **Should be reapplied** at each dressing change - **Level I**

Standards of Care for Peripheral Intravenous Catheters: Evidence-Based Expert Consensus

Judy Thompson, MSNEd, RN, VA-BC

Association for Vascular Access, San Diego, CA

Marlene M. Steinheiser, PhD, RN, CRNI

Infusion Nurses Society, Norwood, MA

J. Blake Hotchkiss, MSN, MBA, RN, VA-BC

MaineHealth Franklin Hospital, Farmington, ME

James Davis, MSN, RN, CCRN-K, HEM, CIC, FAPIC

ECRI, Plymouth Meeting, PA

Michelle DeVries, MPH, FAPIC, CPHQ, CIC, VA-BC

ICU Medical, San Clemente, CA

Katie Frate, BSN, RN, VA-BC

St. Joseph's Children's Hospital, Tampa, FL

Robert Helm, MD

Portsmouth Regional Hospital, Portsmouth, NH

Chris W. Jungkans, BSN, BS, RN, VA-BC, CPUI

Advocate Health, WI

Swapna Kakani, MPH

Swapna Kakani Consulting, LLC, Huntsville, AL

Sean Lau, MBA, BSN, RN, VA-BC

Stanford Health care, Palo Alto, CA

Karen Lindell, DNP, RN

Moore Regional Hospital and FirstHealth of the Carolinas, Pinehurst, NC

Kristen McNiff Landrum, MPH

KM Healthcare Consulting, Wilmington, NC

Karen A. McQuillan, MS, RN, CNS-BC, CCRN, CNRN, TCRN, FAAN

American Association of Critical-Care Nurses, Aliso Viejo, CA

R Adams Cowley Shock Trauma Center, University of Maryland Medical Center, Baltimore, MD

DJ Shannon, MPH, CIC, VA-BC, FAPIC

IU Health Adult Academic Health Center, Indianapolis, IN

Indiana University Richard M. Fairbanks School of Public Health, Indianapolis, IN

Lorelle Wuerz, PhD, MSN, RN, VA-BC, NEA-BC

New York-Presbyterian Hospital, New York, NY

Stephanie Pitts, DNP, RN, CPN, VA-BC, NEA-BC

B. Braun Medical Inc. Bethlehem, PA

Highlights

- The United States purchases an estimated 350 million PIVCs annually.
- PIVC insertion is the most frequently performed invasive procedure in healthcare.

PIVC Standards

Endorsing Organizations



RCT: Effect of Novel Peripheral Intravenous Catheter Securement for Children to Prevent Catheter Failure

POPULATION

189 Males, 194 Females



Children in emergency departments (EDs) anticipated to require both a peripheral intravenous catheter (PIVC) and admission to hospital for >24 hours

Median age, 36 mo

SETTINGS / LOCATIONS




2 regional hospitals in Queensland, Australia

INTERVENTION

402 Patients randomized

383 Patients analyzed



134 Standard care

Application of a bordered polyurethane dressing over the PIVC

118 Integrated securement dressings (ISD)

Application of an ISD with variable sizes over a PIVC

131 ISD with tissue adhesive (TA)


Use of TA at PIVC site and hub prior to integrated securement dressing

PRIMARY OUTCOME

The primary outcome was PIVC failure, defined as unplanned PIVC removal where PIVC reinsertion was required

FINDINGS

PIVC failure was lowest in ISD with TA group compared with ISD and standard care.



43 Of 134 patients had PIVC failure; IRR reference
24 Of 118 patients had PIVC failure; IRR, 0.80; 95% CI, 0.47-1.35
15 Of 131 patients had PIVC failure; IRR, 0.45; 95% CI, 0.23-0.84

Charters B, Foster K, Lawton B, et al. Effect of Novel peripheral intravenous catheter securement for children and catheter failure: a randomized clinical trial. *JAMA Pediatr*. Published online March 25, 2024. doi:10.1001/jamapediatrics.2024.0167


© AMA

EB4: CLABSI Reduction Through Use of a Tissue Adhesive in a Cardiac Intensive Care Unit

*“Through the successful implementation of a TA in the cardiac ICU, an **88% CLABSI reduction** was achieved. Furthermore, the **dressings-change frequency** was reduced by more than 50%, and no patients were reported as having skin-related reactions.”*

- Rodriguez H, Nuila C, Houston Methodist Hospital


Quality and Safety Performance: Where we are Now!




Updated data presented at AVA 2022:
Repeat of results following TA Use Pause due to COVID



PIVC Sealed and Secured with Tissue Adhesive



Tissue Adhesive Securement: Insertion Site Sealed







Proven Clinical Benefits

Reduced insertion site bleeding and oozing: 80%
Reduced dressing changes: 45% to 57%

Reduced migration and dislodgement: 93%

Seals and protects the insertion site to reduce the risk of infection

Reduce CLABSI Risk: 50% to 89%

Reduce PIV Failure Risk: 63%

CVC/PICC COST SAVINGS

50% Reduction in PICC/CVC Dressing Changes = > SecurePortIV® directly pays for itself

Additional PICC/CVC Financial Benefits

93% Reduction in migration

75% Reduction in dislodgement

50% to 89% Reduced risk for CLABSI

Potential to eliminate the use of other engineered securement and antimicrobial dressings

PERIPHERAL IV COST SAVINGS

74% reduction in PIV dressing disruption and 87% IV therapy completion rate

65% Reduction in PIV Failure Risk

SecurePortIV®
Catheter Securement Adhesive
Highly Purified Medical Cyanoacrylate

SecurePortIV® Pays for Itself in Dressing Change Reduction

Dressing change reductions demonstrated in the published literature: **45% to 57%**



Infection Prevention

Reduced risk for
Infection



Reduce Catheter
Replacements

93% Reduction in
catheter migration



Reduce Catheter
Dislodgement

75% Reduction in
catheter
dislodgement



Eliminate Other
Products

Engineered securement
Hemostatic Products
Antimicrobial disk or gel



Learn more
about SecurePortIV®

NQF Hospital-Onset Bacteremia (HOB) Playbook



THANK YOU

Jennifer.Perou@HBFuller.com

562.619.3022

Jennifer Perou, RN
Regional Business Manager CA & HI

- Charters, B., Foster, K., Lawton, B., Lee, L., Byrnes, J., Mihala, G., Cassidy, C., Schults, J., Kleidon, T. M., McCaffery, R., Van, K., Funk, V., & Ullman, A. (2024). Novel Peripheral Intravenous Catheter Securement for Children and Catheter Failure Reduction: A Randomized Clinical Trial. *JAMA Pediatr.* <https://doi.org/10.1001/jamapediatrics.2024.0167>
- H. B. Fuller Medical Adhesive Technologies. SecurePortIV: Indication for Use (IFU). Retrieved August 1, 2024, from www.SPIVTraining.com
- Kleidon, T. M., Ullman, A. J., Gibson, V., Chaseling, B., Schoutrop, J., Mihala, G., & Rickard, C. M. (2017). A Pilot Randomized Controlled Trial of Novel Dressing and Securement Techniques in 101 Pediatric Patients. *J Vasc Interv Radiol*, 28(11), 1548-1556 e1541. <https://doi.org/10.1016/j.jvir.2017.07.012>
- National Quality Forum (NQF). (September 2024). Hospital Onset-Bacteremia and Fungemia Playbook. Retrieved October 1, from https://www.qualityforum.org/Publications/2024/9/Hospital-Onset_Bacteremia_and_Fungemia_Playbook.aspx
- Nickel, B., Gorski, L., Kleidon, T., Kyes, A., DeVries, M., Keogh, S., Meyer, B., Sarver, M. J., Crickman, R., Ong, J., Clare, S., & Hagle, M. E. (2024). Infusion Therapy Standards of Practice, 9th Edition. *Journal of Infusion Nursing*, 47(1S), S1-S285. <https://doi.org/10.1097/nan.0000000000000532>
- Prince, D., Kohan, K., Solanki, Z., Mastej, J., Prince, D., Varughese, R., & Patel, M. (2017). Immobilization and Death of Bacteria by Flora Seal® Microbial Sealant. *International Journal of Pharmaceutical Science Invention*, 6(6), 45-49. [http://www.ijpsi.org/Papers/Vol6\(6\)/G06064549.pdf](http://www.ijpsi.org/Papers/Vol6(6)/G06064549.pdf)
- Ralph Webber, J. L., & Maningo-Salinas, M. J. (2020). “Sticking It to Them”—Reducing Migration of Peripherally Inserted Central Catheters. *Journal of the Association for Vascular Access*, 25(1), 10-15. <https://doi.org/10.2309/j.java.2020.001.004>
- Rodruiguez H., & Nuila C. (2022). National Teaching Institute Evidence-Based Solutions Abstracts. *Crit Care Nurse*, 42(2), e11-e12. <https://doi.org/10.4037/ccn2022820>
- Thompson, J., Steinheiser, M. M., Hotchkiss, J. B., Davis, J., DeVries, M., Frate, K., Helm, R., Jungkans, C. W., Kakani, S., Lau, S., Lindell, K., Landrum, K. M., McQuillan, K. A., Shannon, D. J., Wuerz, L., & Pitts, S. (2024). Standards of Care for Peripheral Intravenous Catheters: Evidence-Based Expert Consensus. *Journal of the Association for Vascular Access*, 29(3), 15-26. <https://doi.org/10.2309/java-d-24-00011>
- Ullman, A. J., Kleidon, T., Gibson, V., McBride, C. A., Mihala, G., Cooke, M., & Rickard, C. M. (2017). Innovative dressing and securement of tunneled central venous access devices in pediatrics: a pilot randomized controlled trial. *BMC Cancer*, 17(1), 595. <https://doi.org/10.1186/s12885-017-3606-9>
- van Rens, M., Nimeri, A. M. A., Spencer, T. R., Hugill, K., Francia, A. L. V., Olukade, T. O., & Mahmah, M. A. (2021). Cyanoacrylate Securement in Neonatal PICC Use: A 4-Year Observational Study. *Adv Neonatal Care*. <https://doi.org/10.1097/ANC.0000000000000963>
- Zhang, S. (2023). *White Paper: Effective Immobilization of Candida auris by SecurePortIV*. HB Fuller | Adhezion Biomedical.



H.B. Fuller

| Medical Adhesive Technologies